

## Stickpack filling and sealing machine

High-precision dosing into Stickpacks



# **MFS800**





### Range of products:

### Liquids

e.g. vaccines, reagents, analysing liquids...

#### **Powders**

e.g. pharmaceutical powders...

### **Granules, small piece goods**

e.g. pharmaceutical products, micro tablets  $\dots$ 









### Stickpack filling and sealing machine

### Type: MFS800





### Description of operation

This system is mainly developed to control the net weight of a Stickpack by leaving out the influence of the weight of the wrapping material.

To this end, empty bags are put into carriers - manually or supplied by a MERZ packaging machine.

The filling and weighing of the bags is performed in several steps. The filled Stickpack carriers are being conveyed step by step via a transport system.

- 1. Step: Additional Stickpack spreading device
- 2. Step: The empty bag is weighed together with the carrier (tara)
- 3. Step: Additional camera for controlling open Stickpacks
- 4. Step: Dosing or rather filling of product 1
- 5. Step: Weighing of product 1
- 6. Step: Horizontal sealing product 1 (Option for separation of product 1 and product 2)
- 7. Step: Dosing or rather filling of product 2
- 8. Step: Weighing of product 2
- 9. Step: Horizontal sealing (closing of Stickpack)
- 10. Step: Waiting position (reserve position e.g. for control by camera)11. Step: Sachet sorting into good and bad bags with reject verification
- 12. Step: Sensor with holder after Stickpack transfer for check of empty workpiece carriers
- 13. Step: The empty carrier moves on, back to the starting position

### **Technical Details**

- Machine completely out of stainless steel and aluminum
- Pneumatic parallel sealing drive
- Sealing temperature left/right separately adjustable
- Bag counter
- Improved control by determination of empty bag net weight
- Siemens PLC

#### **Options**

- Discharge conveyor
- Dust collection
- Peel seam between two chambers possible
- Possibility of filling of two different products into one bag
- Transformer for special voltage

#### **Technical Data**

Bag width	9-23 mm
Bag length	30-120 mm
Filling volume*	ranging from some microliters
	up to 5ml or up to 200 micro tablets
Output bags/min	up to 25
Air connection	6 bar
Wattage	approx. 2 KW
Voltage supply standard	400V/3/N/PE/50Hz
Protection class	IP 54
Weight	approx. 1150 kg
Colour standard	stainless steel - alu

\*Depending on product, filling quantity and packaging material



